



LOW RETENTION FILTERED PIPETTE TIP

Part Numbers:

10 uL Pipette Tip | MLS-PTFS-10L
96 tips per rack, 960 tips per pack

20 uL Pipette Tip | MLS-PTFS-20L
96 tips per rack, 960 tips per pack

200 uL Pipette Tip | MLS-PTFS-200L
96 tips per rack, 960 tips per pack

1000 uL Pipette Tip | MLS-PTFS-1000L
96 tips per rack, 768 tips per pack

Manufacturer: LabCon

Color: Natural

Sterile: Yes

Filtered: Yes

Fit: Universal

Material: Polypropylene

Certified free: RNase/DNase, ATP, endotoxin/pyrogens (<0.05EU/mL), cytotoxic effects, and trace heavy metals (<1ppm).

