

MLS-GL-981

Blue Nitrile Gloves



Description

- Examination gloves
- Ambidextrous
- Beaded Cuff - providing additional strength and preventing liquid roll off
- LatexFree- to reduce the risk of hypersensitivity and allergic reactions
- Nitrile Soft Gloves are a softer and thinner for improved feel and sensitivity, with Micro Textured fingers for excellent wet and dry grip
- Superior barrier protection over latex and vinyl
- Excellent puncture resistance and strength
- Powder Free
- 100 gloves per pack, 10 packs per case

Dimensions

Palmsize	Palm Width
S	80 mm
M	95 mm
L	110 mm
XL	110 mm

Test results

Test type	Inspection level	AQL
Watertight test	G1	1.5
Watertight test	G1	2.5
Minor defects	G1	4.0
Dimensions	S2	4.0
Physical Properties	S2	4.0

Thickness

Measurement type	Single wall (mil)
Finger	Min 4.0
Palm	Min 3.1
Cuff	Min. 2.4

Keyfacts

Accreditations:

Medical Device Directive -93/42/EEC Class 1 (2007/47/EC)

EN455-1,2,3,4 - Medical Gloves for Single use.

EN420; Protective Gloves; General Requirements

EN374-1,2,3, 4 & 5 Determination of Resistance to Water Penetration.

Leak test passed.

ISO 9001
ISO 13485

Conforms to ASTM D6978

Physical Properties

Property	Before Aging	After Aging
Tensile strength	Min.14	Min.14
Elongation (%)	Min.650	Min.500
Force at Break	6.0nm	6.0nm

Symbols:



Approved by:

